



# LightSheer® DUET™

High Speed Permanent Hair Reduction

RETHINK THE BUSINESS OF HAIR REMOVAL



HIGH SPEED


# Introducing LightSheer® DUET™

Rethink the Business of Hair Removal


We've revolutionized the Gold Standard LightSheer with vacuum-assist technology and by dramatically increasing the spot size to 22 x 35 mm. The LightSheer DUET enables you to **treat backs and legs comfortably** in only 15 minutes, making hair removal profitable...**without anesthetics** and gels.

- No topical anesthetic required
- Eliminate prep time
- Reduce treatment times up to 75%
- Differentiate your practice
- Grow your patient base
- Grow your profits


## High-speed integrated vacuum technology



1) vacuum draws skin into hand-piece; 2) dermis is stretched thin; 3) hair follicle is pulled closer to energy source



4) melanocytes are spread apart; blood flow is constricted; 5) energy is applied




6) skin is released; 7) hair follicles are damaged; 8) tip moved to next location

- Reduction in competing chromophores enables more efficient delivery of laser energy
- Lower required energy levels produce a more comfortable procedure, eliminating the need for topical anesthetics

## depth of penetration

Larger spot size enables greater depth of penetration, for improved absorption of energy at the target follicle<sup>1</sup>



"I've used several laser hair removal devices in my practice... LightSheer remains the 'gold standard' because it is the most effective for targeting unwanted hair. Effectiveness, long lasting results, satisfied patients...all add up to repeat business and positive referrals."

**Beth Haney, NP**

"Patients will love it. Nurses will love it... I recently treated a large gentleman's back and in the 20 minutes it took to treat him, he was practically falling asleep – the treatment was that comfortable."

**Christine Dierickx, MD**



# er

## Versatile

The LightSheer DUET includes two **optimal wavelength 800 nm** diode hair removal systems on one versatile platform.



### LightSheer ET

9 x 9 mm spot size

- Fluence up to 100 J/cm<sup>2</sup>
- ChillTip™ continuous contact cooling
- Pulse widths up to 400 ms


NEW

### LightSheer HS

22 x 35 mm spot size

- Largest spot size of any diode
- Integrated vacuum-assist technology
- Pulse widths up to 400 ms
- No cooling or anesthesia required

actual spot size



Intuitive touchscreens

- Recommended treatment settings
- Customizable treatment presets

## Differentiate your practice

- Treat with 1/10 the pulses—save time and effort
- Make large-area treatments affordable for your patients

Estimated Pulses per Treated Site

Treatment Site	Est. Area (cm <sup>2</sup> )	Spot Size:		
		9x9mm	12x12mm	22x35mm
bikini area	125	154	87	16
underarm (2)	190	235	132	25
man's back	2750	3395	1910	357
full legs (2)	8100	10,000	5695	1052

\*12 x 12 mm spot size available on XC platform only

"The new LightSheer DUET is extremely well tolerated by patients. My hair removal nurse likes the single pulsing for treatment speed and patient comfort. My expectation is the LightSheer diode will perform better as it penetrates deeper...This system is definitely worth pursuing."

**Girish Munavalli, MD**

"I had my back treated with the new High Speed system in 15 minutes. Normally, this would take about an hour to do, so the time savings alone is significant. ...I think it's a winner."

**Mitch Goldman, MD**

## configurations

- **LightSheer Duet with HS** – Exciting new alternative for those seeking high-speed, comfortable treatments to offer your patient population, and increase your patient throughput.
- **LightSheer Duet with ET & HS** – Versatile platform for those focused on practice growth, enables a customized approach across the broadest range of patients and treatment areas.

## specifications

PLATFORM	LightSheer DUET	
SYSTEMS AVAILABLE	LightSheer ET	LightSheer HS
<b>Technical Specifications</b>		
Laser Type	AlGaAs Diode Arrays	AlGaAs Diode Arrays
Wavelength	800 nm (nominal)	800 nm (nominal)
Spot Size	9 mm x 9 mm	22 mm x 35 mm
Repetition rate	Up to 3 Hz	Up to 3 Hz
Pulse Repetitions	1 pulse	1 – 3 pulses
Fluence	10 - 100 J/cm <sup>2</sup>	4.5 - 12 J/cm <sup>2</sup>
Pulse Width (OptiPulse™)	5 - 400 ms	30 - 400 ms
Peak Power	1600 W	1600 W
Comfort Management	ChillTip™ contact cooling (5°C)	Vacuum-Assist (up to 18 inHg)
Preset Parameters	Yes	Planned
Warranty	1 Year	1 Year
On-Site Service	Yes; replaceable handpiece	Yes; replaceable handpiece

### Platform Specifications

Electrical	100-240 VAC, 15 A max. 50/60 Hz single-phase grounded outlet
Weight	48 kg (105 lbs)
Dimensions (WxDxH)	44 x 50 x 112 cm; (17.4 x 19.7 x 44 in)

<sup>1</sup> Klavuhn KG. Illumination Geometry: The Importance of Laser Beam Spatial Characteristics. *Lumenis*, 2001.

FDA cleared for hair removal and permanent reduction in all skin types (Fitzpatrick I-VI) and tanned skin; PFB (pseudofolliculitis barbae); superficial leg veins; benign pigmented lesions.



#### AMERICAS

**Lumenis Inc.**  
Santa Clara, CA, USA  
Tel +1 408 764 3000  
1 877 586 3647  
Fax +1 408 764 3999  
1 800 505 1133

#### EUROPE / MIDDLE EAST

**Lumenis (GERMANY) GmbH**  
Dreieich-Dreieichenhain,  
Germany  
Tel +49 6103 8335 0  
Fax +49 6103 8335 300

#### JAPAN

**Lumenis (JAPAN) Co. Ltd.**  
Tokyo, Japan  
Tel +81 3 5789 8300  
Fax +81 3 5789 8310

#### ASIA / PACIFIC

**Lumenis (CHINA) Ltd.**  
Beijing, China  
Tel +86 10 6510 2620  
Fax +86 10 6510 2621